




dragonfli

koppert.com



Koppert

Spidex
Vital Plus



Spidex Vital Plus (*Phytoseiulus persimilis*) - Premium Sachet Application Instructions

Please use immediately upon receipt. For indoor use only.

Required conditions:

- Temperature range: between 13 and 27°C
- Humidity: 70% (lower humidity will not kill the *phytoseiulus* predators but will make reproduction more difficult)
- Do not apply the *phytoseiulus* if insecticides were used in the past 1-2 weeks

Before application:

- If possible, manually remove any webbing on the leaves of your plants
- Prune any heavily damaged leaves
- Do not open the sachets. The sachets already contain an emergence hole

Application:

- Use all of the sachets at once and do not store
- Hang the sachets onto stems, or leaves, at a rate of 1-2 per plant
- Hang the sachets out of direct sunlight
- Do not use any insecticides once the predators have been applied
- The sachets will release predators for 2 weeks

Next Steps:

- Observe Spider Mite activity during the 2 week life cycle of the sachet and ensure to replace the sachets to continue protection
- The sachets can be used with our ***Amblyseius californicus*** predators
- If Spider Mite numbers rapidly increase switch to a curative application of our ***Phytoseiulus persimilis***

www.dragonfli.co.uk